
haze-gard i Pro

Short Instructions

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Instructions rapides

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Instrucciones resumidas

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Table of Contents

English 4

Deutsch..... 14

Français..... 24

Italiano 34

Español 44

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
Table of Contents

- 1 System Description 5**
 - 1.1 Preparation..... 5
 - 1.2 Connectors..... 5
 - 1.3 Operation 6
 - 1.4 Cleaning..... 6
- 2 Software Installation 7**
 - 2.1 Installation..... 7
 - 2.2 Configuration..... 7
- 3 Display and Main Menu 8**
- 4 Quick Measurements 9**
- 5 Standard Management 10**
- 6 Project Management..... 11**
- 7 System Configuration 12**
- 8 Technical Data 13**

1 System Description

The entire system consists of instrument, black cover, power supply, foot switch and software. The instrument can be operated in horizontal or vertical orientation.


1.1 Preparation

1. Connect power supply and foot switch to instrument.
2. Connect instrument with PC the USB cable type A - B.
3. Download and install "smart-chart" software.
4. Remove protective covers from clarity and haze port.
5. Turn instrument on with power switch.
6. Run initial configuration wizard: Select language, enter date / time and choose a measurement method (**ASTM Illuminant A**, **ASTM Illuminant C** and or **ISO**).
7. A calibration is necessary at first startup. You will be guided by the instrument. For manual calibration - recommended after warm-up before collecting measurements - select **Configuration > Calibration**.

1.2 Connectors

Connectors and power switch are placed on the rear side.

- | | |
|-----------------------------|------------------------------|
| 1 Connector for foot switch | 2 LAN interface |
| 3 USB type A interface | 4 Connector for power supply |
| 5 Power switch | |


NOTICE

Retain the packaging in case the instrument needs to be shipped.

1.3 Operation

To ensure a repeatability and reproducibility of 10% of measured value but > 0.01 units in a haze range $\leq 0.3\%$, the following conditions must be strictly applied:

- Warm up the instrument for 60 minutes prior to usage.
Note: It is recommended to keep the instrument switched on over night.
- Perform a calibration after warm-up time or if used continuously once a day.
- Attach the black, magnetic cover to ensure the measurement port is in the shade – a horizontal instrument set-up is recommended.

- Observe standard operating conditions for electronic instruments:
 - Avoid excessive oscillations and vibrations.
 - Avoid extreme ambient temperatures or rapid changes in temperature.
 - Avoid relative humidity more than 85%.
 - Avoid extreme electromagnetic fields.

1.4 Cleaning

To verify the accuracy of the measurement results a checking standard close to the product's haze level is recommended. For very low haze applications the haze standard $< 0.3\%$ (4798) is recommended.

- A regular verification, once a week, is advisable.
- If the measured value of the checking standard is not within the defined range printed on the standard's label, proceed as follows:
 - Illuminate with a strong flashlight and tilt the standard approx. 45° to see smallest dust particles:

- For cleaning use a fresh, one way, lint-free cloth with some isopropanol. Wipe only the glass and not the metal frame in a circular movement with low pressure.

2 Software Installation

2.1 Installation

1. Download zip-file from:
<https://www.byk-instruments.com/haze-gardi>
2. Save the file into a new folder and extract the complete archive.
3. Right mouse click on "**install.exe**" and select option "**Run as administrator**".
4. Follow the setup instructions on the screen.


After installation the software package "smart-lab" can be used for **30 days** free trial. Thereafter, you need to register the package.

2.2 Configuration


In smart-lab go to **Configuration > Options**.

- Set **Language**.
- Define **Autosave** interval (option).
- Activate **Precision mode haze-gard i** to display haze results with 3 digits.

3 Display and Main Menu

The display unit provides following elements.

1	Touch screen Tap the icons on the screen to in- voke the various functions.	2	Status LED Indicates the current status during measurement.
3	USB interface Connect USB flash drive to transfer data, update firmware or save ser- vice file.	4	Measure button Trigger a measurement, see " Quick Measurement [9]".


The main menu provides following elements.

 Measure [9] Enter quick measurement mode.	 Projects [11] Manage projects with test series and save measurement data automatic- ally.
 Standards [10] Define test method (scales) and tol- erance values for pass / fail indica- tion.	 Configuration Check calibration and change instru- ment settings.
 System Time Displays the current time, see Con- figuration > Set Date / Time.	 Shut Down To avoid loss of data or malfunctions always tap the power off icon be- fore switching off the instrument.

4 Quick Measurements

Perform quick measures on samples. Measurements can be triggered via the foot switch, the button **Measure** or the symbol **Measure** on the touch screen. Results can be saved optionally in projects.

1. In the main menu select **Measure**.
 2. Select test method on top; corresponding scales are displayed.
 3. Place sample at appropriate port; depending on scales given in standard:
 - Left port: Measure **Clarity (C)**.
 - Right port: Measure **Haze (H)** and **Total Transmission (T)**.
 4. Press the **Measure** button.
 5. Continue with next measurement or proceed to **Next Sample**.
 6. To store results select **Save**; systems jumps over to a "Project [[11](#)]".
 7. To end quick measurement mode, select the **Home** icon.
- In top row the average value is shown, single readings below.


NOTICE

In case of statistics evaluation the average value will be saved. Number of measurements per sample depends on statistics settings in "[Standards \[10\]](#)" and in "Configuration".

5 Standard Management

Define test method (scales) and tolerances for pass / fail indication (values). You can change existing standards and create new standards.

To create a new standard:

1. In the main menu select **Standards**.
2. In the screen **Standards** select **+New**.
3. Input a **Name** and select the **Scales** required in the standard.
4. In upper right corner select number of measurements for statistics evaluation, e.g. "n = 3".
5. Select **Edit** to change the tolerance value(s) for a scale.
6. To store the standard in memory select **Save**.
7. To end standard management, select the **Home** icon.

To view or change an existing standard: Select **(Various Standards)** and select the standard from the list.

NOTICE

Standards created in "smart-chart" (stored in a separate group) are protected. These standards cannot be edited.

6 Project Management

Open an existing project or create a new one, measure samples and store data automatically. A project consists of standard, test series and samples; path to a measurement is given in top row. Tap on a part of the path to go to that level.

The screenshot displays the software interface for project management. At the top, a breadcrumb trail shows the navigation path: **Foils** > **Product BB8** > **Batch 2206** > **SAMPLE 008**. The time is 10:42 and the page is 3/3. Below this, the measurement data is organized into three columns: **ASTM - C**, **H**, and **BYK C**. The current values are **81.9**, **11.4**, and **87.3** respectively. A tabbed interface shows **VALUES** and **STATISTIC**. The **VALUES** tab is active, displaying a table of readings. The table has three columns for different measurement types and a date/time column. The data is as follows:

READING	ASTM - C	H	BYK C	Date/Time
READING 01	78.5	13.2	87.1	2022.08.05 10:41
READING 02	80.8	8.93	86.2	2022.08.05 10:41
READING 03	86.3	12.0	88.5	2022.08.05 10:41

At the bottom, there are buttons for **DELETE**, **RENAME**, and **NEXT SAMPLE**.

1. In the main menu select **Projects**.
2. In the screen **Projects** select **+New**.
3. Input a **Name** for the new project.
4. Load an existing standard or create a new one; edit the standard if needed.
5. Select an existing test series or create a new one.
6. Add measurements to your new project.
7. To end project management, select the **Home** icon.


NOTICE

The color for each measurement indicates the **Pass / Fail** result according to the value(s) entered in the standard.

7 System Configuration

The configuration menu provides following elements.

Calibrate Instrument

Starts the procedure for the system calibration.

Set Date / Time

Set date, time and timezone.

Continuous Measurements

The "Continuous Measurement" symbol appears in the main menu.

Single Measurement Mode

Shortened measurement time for automation; can only be used in a dark environment.

Select Language

Select preferred display language.

Rotate Display

Set the display to horizontal or vertical orientation.

Themes

Switch the color scheme on the display (light/dark).

LAN Settings

Check and change network configuration.

System Info (*)

View Serial No. / Catalog No. / Firmware Version.

Statistics

On: Number of measurements according to standard is made and evaluated.

Off: One single measurement per sample is made.

Statistics Fixed

On: Exact number of measurements according to standard is to be made.

Off: A free number of measurements can be made.

(*) Factory reset can be performed under **System Info**. There is an option to delete all standards / projects. Save these data to a USB flash drive before.

8 Technical Data

Instrument

Meas. geometry	0° / Diffuse
Sample port	25.4 mm (0.85 in)
Measurement area	18.0 mm (0.65 in)
Spectral response	CIE luminosity function y
Illuminants	ASTM-D1003, CIE-A, CIE-C, CIE-D65, ISO 13468, 14782
Measurement range	Haze: 0 - 100% Clarity: 0 - 100% Total transmission: 0 - 100%
Repeatability (standard deviation)	± 0.1 units For Haze ASTM-C < 0.3 %: ± 10% of measured value or at least > 0.01 units
Reproducibility (standard deviation)	± 0.4 units For Haze ASTM-C < 0.3 %: ± 10% of measured value or at least > 0.01 units
Statistics memory	5000 readings
Interfaces	LAN USB 2.0 Additional front USB port for flash drive

General

Temperature range	10° C to 40° C (50° F to 104° F) for operation 0° C to 50° C (32° F to 122° F) for storage
Relative humidity	Up to 85% non-condensing at 35° C (95° F)
Operation altitude	Up to 2000 m (6561 ft)
Dimensions (LxWxH)	62 x 33 x 22 cm (24 x 13 x 9 in)
Weight	18 kg (40 lbs)
Languages	English, German, French, Italian, Polish, Spanish, Chinese, Japanese, Russian
Power supply	Input: 12 V  , max. 600 mA
External power supply	Input: 100 - 240 V  , 50 - 60 Hz, 620 mA Output: 12 V  , max. 2.5 A Vendor: Adaptertech Model: ATS030 Series Safety approvals: CB / UL / cUL / FCC / GS / CE / PSE / BSMI / CCC / RCM / LPS

Technical data are subject to change without notice.

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- 1 Systembeschreibung..... 15**
 - 1.1 Vorbereitung 15
 - 1.2 Anschlüsse..... 15
 - 1.3 Bedienung 16
 - 1.4 Reinigung 16
- 2 Softwareinstallation..... 17**
 - 2.1 Installation..... 17
 - 2.2 Konfiguration..... 17
- 3 Anzeige and Hauptmenü 18**
- 4 Schnelle Messungen 19**
- 5 Standardmanagement..... 20**
- 6 Projektmanagement 21**
- 7 Systemkonfiguration 22**
- 8 Technische Daten 23**

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BYK-Gardner GmbH