byko-visc RT Lite

Short Instructions

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Instructions rapides

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Instrucciones resumidas

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Table of Contents

English 4

Deutsch..... 14

Français..... 24

Italiano 34

Español 44

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Table of Contents

1	System Description	5
2	Instrument Setup	6
3	Main Menu	7
4	Measurement Settings	8
5	System Configuration	9
6	Viscosity Measurement	10
7	Browse Measurements	11
8	Technical Data	12
9	Service Points	13

1 System Description


The entire system consists of instrument, spindle and power supply.

CAUTION!


Read the safety instructions and the full instruction manual - available via download - before using this instrument. The viscometer has moving parts which should never be touched during operation.

Overview of Parts

Viscometer Head


NOTICE

Always remove the spindle and attach the protective cap before moving the instrument. Retain the packaging in case the instrument needs to be shipped.

2 Instrument Setup


1. Place the viscometer in a suitable location.
2. Turn protective cap clockwise to access spindle attachment.
3. Plug the power supply into a mains outlet and plug the output cable into the back of the viscometer head.
4. Turn instrument on with the power button located on the back of the viscometer head.
5. After power on a level check takes place. Use the adjustable feet to level the instrument.
6. A message appears when instrument is leveled.

After finishing the leveling the instrument will automatically go to the main menu.

3 Main Menu

In the main menu the following symbols are displayed.

1 Viscosity Measurement

Perform measurements using existing measurement settings.

2 Browse Measurements

Manage measurement results of last measurement performed.


3 Measurement Settings

Change measurement settings and / or instrument configuration.

4 Measurement Settings

Tap the settings icon in the main menu to manage settings changed **often**.

Spindle	Select spindle based on spring. Only spindles applicable to your model will be shown.
Speed	Select a speed from the predefined list from 0.3 to 200 rpm (in 21 steps). Spindle and speed combination will show the maximum viscosity that can be measured with this combination.
Density	Available if Density is activated in System Configuration .
Time to Stop	Set measurement time in hours / minutes / seconds.
Time to Torque	Set torque measurement at which the unit will stop when reached.
Configuration	Check and change system settings, see System Configuration .

Use back navigation icon (<) to go back to previous page. Current settings will be saved.

NOTICE

Settings for **Speed** and **Spindle** will be used as defaults and can also be changed directly during measurements.

5 System Configuration

Tap on the configuration icon in the **Measurement Settings** menu to manage settings changed **seldom**.

<

Configuration

 Units

mPa.s, °C

 Density

☐

 Language

English

 Date / Time

17 17

 Color Scheme

☐

 Information

 Factory reset

 Zero Check

 Lock Settings

☐

Units	Select units for display of viscosity and temperature.
Density	Turn option Density on or off in Measurement Settings .
Language	Select preferred display language.
Date / Time	Set date, time and display format. Will be used for the Browse function.
Color Scheme	Change color theme on display (bright / dark).
Information	View current firmware version and serial number for instrument.
Factory Reset	Fall back to factory settings.
Zero Check	Perform a zero check. For details see operating manual.
Lock Settings	Lock current settings against unintentional changes.


NOTICE

Configuration for **Units** will be used as default and can also be changed directly during measurements.

6 Viscosity Measurement

1. If protective cap is attached, turn it clockwise to remove.
2. Move spindle attachment upwards, insert spindle and release attachment.
3. Tap measurement icon in main menu to access measurement screen.
4. Lower unit into sample and tap green **Start** button.


1	TTS / TTT	Counter for Time to Stop or Time to Torque
2	Viscosity	Average of measurements - taken every spindle revolution, tap unit to change
3	Maximum	Maximum value, depending on current spindle and speed
4	Speed	Current spindle speed in rotations per minute (RPM), tap value to change
5	Spindle	Identifier for current spindle, tap entry to change
6	Torque	Current value in percent: Values < 15% or > 95% indicate that speed or spindle settings should be changed (invalid measurement)
7	Temp.	Current temperature of sample in °C or °F (if temperature probe is attached), tap to change unit
8	Graph	Live updating graph of viscosity over time
9	Start	Green Start button: Start and resume measurement
10	Pause	Blue Pause button: Pause measurement
11	Stop	Red Stop button: End measurement and save results
12	Help	White Question mark: Open context-sensitive help

Touch areas of the screen with measured data to access corresponding setting.

7 Browse Measurements

Tap the browse icon in the main menu to manage saved measurement data. The final viscosity value for a measurement is displayed here. Export to USB to see all viscosity values measured over time.

1	File Name	Identify a test series by its name, tap to rename.
2	Final Visc	Gives the final viscosity value for this test series.
6	Select All	Select / deselect all rows / test series at once.
3	Write USB	Write selected test series to USB flash drive.
4	Print USB	Print selected test series to USB printer.
5	Delete	Delete selected test series from memory.

Active options are dark blue. Inactive options are light blue.

8 Technical Data

Instrument

Catalog number	8350 byko-visc RT Lite L (Low viscosity) 8351 byko-visc RT Lite R (Regular viscosity) 8352 byko-visc RT Lite H (High viscosity)
Measurement range	byko-visc RT Lite L: 1* - 2,000,000 cP byko-visc RT Lite R: 5* - 13,000,000 cP byko-visc RT Lite H: 300 - 106,000,000 cP * With Low Viscosity Adapter (LVA) accessory
Spindles included	byko-visc RT Lite L: 4 byko-visc RT Lite R: 6 byko-visc RT Lite H: 6
Speed range	0.3 to 200 rpm in 21 steps
No. of speeds	21
Repeatability	0.2%
Accuracy	1.0% of range
Measurement features	Density, Single-point Tests
Programmable features	Time to Stop, Time to Torque
Temperature measurement	PT100 temperature probe (optional)
Display	Capacitive color display
Languages	English, German, French, Italian, Spanish, Russian, Japanese, Chinese

General

Temperature range	10° C to 40° C (50° F to 104° F) for operation 0° C to 60° C (32° F to 140° F) for storage
Relative humidity	Up to 85% non-condensing at 35° C (95° F)
Operation altitude	Up to 2000 m (6561 ft)
Dimensions (LxWxH)	250 x 310 x 470 mm (9.8 x 12.2 x 18.5 in)
Interface	2x USB Type A
Weight	6.5 kg (14.3 lbs)
Power supply	24 V  , max. 1.25 A
External power supply	Input: 100 - 240 V  ; 50 - 60 Hz ; max. 0.7 A Output: 24 V  ; max. 1.25 A Vendor: Mean Well Model: GE30I24 Safety standards: UL62368-1, CSA C22.2 NO. 62368-1, TUV BS EN/EN62368-1, AS/NZS 60950.1, CCC GB4943, EAC TP TC 004 approved

Technical data are subject to change without notice.

9 Service Points

Global service centers with ISO / IEC 17025 accredited laboratories

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Service Point Austria, Hungary, Slovenia

c/o Friedrich W. Bloch
GmbH
Wagramerstrasse 201, 1210
Vienna, Austria

Service Point India

c/o IMCD India
Private Limited
1101-03, B-Wing, ONE BKC,
Bandra Kurla Complex,
Bandra East, Mumbai, MH.
Pin.: 400 051. India

Service Point China

c/o BYK (Tongling) Co. Ltd.
Shanghai Branch
Block 6A, Building A, No 88
Hong Cao Road, Xuhui Dis-
trict, Shanghai 200233, P.R.
China

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BYK-Gardner GmbH