## byko-visc RT

Short Instructions

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Instructions rapides

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Instrucciones resumidas

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# Table of Contents

English ..... 4

Deutsch..... 14

Français..... 24

Italiano ..... 34

Español ..... 44

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# Table of Contents

|          |                                    |           |
|----------|------------------------------------|-----------|
| <b>1</b> | <b>System Description .....</b>    | <b>5</b>  |
| <b>2</b> | <b>Instrument Setup .....</b>      | <b>6</b>  |
| <b>3</b> | <b>Main Menu.....</b>              | <b>7</b>  |
| <b>4</b> | <b>Measurement Settings .....</b>  | <b>8</b>  |
| <b>5</b> | <b>System Configuration.....</b>   | <b>9</b>  |
| <b>6</b> | <b>Viscosity Measurement .....</b> | <b>10</b> |
| <b>7</b> | <b>Browse Measurements .....</b>   | <b>12</b> |
| <b>8</b> | <b>Technical Data .....</b>        | <b>13</b> |

# 1 System Description


The entire system consists of instrument, spindle and power supply.

## CAUTION!


Read the safety instructions and the full instruction manual - available via download - before using this instrument. The viscometer has moving parts which should never be touched during operation.

## Overview of Parts

## Viscometer Head


## NOTICE

Always remove the spindle and attach the protective cap before moving the instrument. Retain the packaging in case the instrument needs to be shipped.

## 2 Instrument Setup


1. Place the viscometer in a suitable location.
2. Turn protective cap clockwise to access spindle attachment.
3. Plug the power supply into a mains outlet and plug the output cable into the back of the viscometer head.
4. Turn instrument on with the power button located on the back of the viscometer head.
5. After power on a level check takes place. Use the adjustable feet to level the instrument.
6. A message appears when instrument is leveled.

After finishing the leveling the instrument will automatically go to the main menu.

### 3 Main Menu

In the main menu the following symbols are displayed.

**1 Viscosity Measurement**

Perform measurements using existing measurement settings and save results in instrument memory.

**2 Browse Measurements**

View and manage measurement results saved in instrument memory.

**3 Measurement Settings**

Change measurement settings and / or instrument configuration.

## 4 Measurement Settings

Tap the settings icon in the main menu to manage settings changed **often**.

<

Measurement Settings

Spindle

TR-8/SC4-21

Speed

5 RPM

Time to Stop

☐

Time to Torque

☐

Ramp

☐

Limits

☐

Curve fitting

☐

Configuration

1

2

3

4

5

6

|                   |                                                                                                                                          |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Spindle           | Select spindle based on spring. Only spindles applicable to your model will be shown.                                                    |
| Speed             | Input a value from 0.1 - 250 rpm. Cone and speed combination will show the maximum viscosity that can be measured with this combination. |
| Density           | Available if <b>Density</b> is activated in <b>System Configuration</b> .                                                                |
| Time to Stop      | Set measurement time in hours / minutes / seconds.                                                                                       |
| Time to Torque    | Set torque measurement at which the unit will stop when reached.                                                                         |
| Ramp              | Define measurement ramp: Enter start speed, end speed, time to run.                                                                      |
| Limits            | Select alarm criteria: Set upper limit, set lower limit, select alarm type (buzzer, display flash, both, or none).                       |
| Curve Fitting     | Choose between: Newton, Bingham, Casson, Ostwald, Hershel-Buckley to overlay onto the measurement graph.                                 |
| Configuration     | Check and change system settings, see <b>System Configuration</b> .                                                                      |
| Presets (1 ... 6) | Settings are automatically stored in selected slot. Tap a slot to recall stored settings.                                                |

Use back navigation icon (<) to go back to previous page. Current settings will be saved.

**NOTICE**

Settings for **Speed** and **Spindle** will be used as defaults and can also be changed directly during measurements.

## 5 System Configuration

Tap on the configuration icon in the **Measurement Settings** menu to manage settings changed **seldom**.

<

Configuration

Units

mPa.s, °C

Density

☐

Language

English

Date / Time

17:15

Color Scheme

☐

Information

Factory reset

Save

☐

Lock Settings

☐

Zero Check

Custom spindle Entry

|                             |                                                                                                                     |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Units</b>                | Select units for display of viscosity and temperature.                                                              |
| <b>Density</b>              | Turn option <b>Density</b> on or off in <b>Measurement Settings</b> .                                               |
| <b>Language</b>             | Select preferred display language.                                                                                  |
| <b>Date / Time</b>          | Set date, time and display format. Will be used for the <b>Browse</b> function.                                     |
| <b>Color Scheme</b>         | Change color theme on display (bright / dark).                                                                      |
| <b>Information</b>          | View current firmware version and serial number for instrument.                                                     |
| <b>Factory Reset</b>        | Fall back to factory settings.                                                                                      |
| <b>Zero Check</b>           | Perform a zero check. For details see operating manual.                                                             |
| <b>Save</b>                 | Turn saving of results on or off in <b>Viscosity Measurement</b> .                                                  |
| <b>Lock Settings</b>        | Lock current settings against unintentional changes.                                                                |
| <b>Custom Spindle Entry</b> | Manage custom spindles: View list of existing custom spindles. For new spindle enter coefficient and a unique name. |


**NOTICE**  
Configuration for **Units** will be used as default and can also be changed directly during measurements.

## 6 Viscosity Measurement

1. If protective cap is attached, turn it clockwise to remove.
2. Move spindle attachment upwards, insert spindle and release attachment.
3. Tap measurement icon in main menu to access measurement screen.
4. Lower unit into sample and tap green **Start** button.


|    |           |                                                                                                                                 |
|----|-----------|---------------------------------------------------------------------------------------------------------------------------------|
| 1  | TTS / TTT | Counter for Time to Stop or Time to Torque                                                                                      |
| 2  | Viscosity | Average of measurements - taken every spindle revolution, tap unit to change                                                    |
| 3  | Maximum   | Maximum value, depending on current spindle and speed                                                                           |
| 4  | Speed     | Current spindle speed in rotations per minute (RPM), tap value to change                                                        |
| 5  | Spindle   | Identifier for current spindle, tap entry to change                                                                             |
| 6  | Torque    | Current value in percent: Values < 15% or > 95% indicate that speed or spindle settings should be changed (invalid measurement) |
| 7  | Temp.     | Current temperature of sample in °C or °F (if temperature probe is attached), tap to change unit                                |
| 8  | Graph     | Live updating graph of viscosity over time                                                                                      |
| 9  | Start     | Green <b>Start</b> button: Start and resume measurement                                                                         |
| 10 | Pause     | Blue <b>Pause</b> button: Pause measurement                                                                                     |
| 11 | Stop      | Red <b>Stop</b> button: End measurement and save results                                                                        |
| 12 | Help      | White <b>Question</b> mark: Open context-sensitive help                                                                         |

Touch areas of the screen with measured data to access corresponding setting.

# Graph and Table Screen

Swipe left and right to access the other screens.


- 13    Graph        Measurement data in graphical form
- 14    Table        Measurement data in table form

## 7 Browse Measurements

Tap the browse icon in the main menu to manage saved measurement data. The final value for a measurement is displayed here. Export to USB to see all viscosity values measured over time.

|   |                   |                                                       |
|---|-------------------|-------------------------------------------------------|
| 1 | <b>File Name</b>  | Identify a test series by its name, tap to rename.    |
| 2 | <b>Final Visc</b> | Gives the final viscosity value for this test series. |
| 6 | <b>Select All</b> | Select / deselect all rows / test series at once.     |
| 3 | <b>Write USB</b>  | Write selected test series to USB flash drive.        |
| 4 | <b>Print USB</b>  | Print selected test series to USB printer.            |
| 5 | <b>Delete</b>     | Delete selected test series from memory.              |

Active options are dark blue. Inactive options are light blue.

# 8 Technical Data

## Instrument

|                         |                                                                                                                                                                 |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog number          | 8356: byko-visc RT L (Low viscosity)<br>8357: byko-visc RT R (Regular viscosity)<br>8358: byko-visc RT H (High viscosity)                                       |
| Measurement range       | byko-visc RT L: 1* - 6,000,000 cP<br>byko-visc RT R: 5* - 40,000,000 cP<br>byko-visc RT H: 300 - 320,000,000 cP<br>* With Low Viscosity Adapter (LVA) accessory |
| Spindles included       | byko-visc RT L: 4<br>byko-visc RT R: 6<br>byko-visc RT H: 6                                                                                                     |
| Speed range             | 0.1 - 250 rpm                                                                                                                                                   |
| No. of speeds           | Freely selectable within speed range                                                                                                                            |
| Repeatability           | 0.2%                                                                                                                                                            |
| Accuracy                | 1.0% of range                                                                                                                                                   |
| Preset memory           | 6 memory locations for measurement settings presets                                                                                                             |
| Measurement features    | Curve Fitting, Density, QC Limits, Multi-point Tests                                                                                                            |
| Programmable features   | Time to Stop, Time to Torque, Multi-step Ramp                                                                                                                   |
| Temperature measurement | PT100 temperature probe                                                                                                                                         |
| Display                 | Capacitive color display                                                                                                                                        |
| Languages               | English, German, French, Italian, Spanish, Russian, Japanese, Chinese                                                                                           |

## General

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Temperature range     | 10° C to 40° C (50° F to 104° F) for operation<br>0° C to 60° C (32° F to 140° F) for storage                                                                                                                                                                                                                                                                                                                            |
| Relative humidity     | Up to 85% non-condensing at 35° C (95° F)                                                                                                                                                                                                                                                                                                                                                                                |
| Operation altitude    | Up to 2000 m (6561 ft)                                                                                                                                                                                                                                                                                                                                                                                                   |
| Dimensions (LxWxH)    | 250 x 310 x 470 mm (9.8 x 12.2 x 18.5 in)                                                                                                                                                                                                                                                                                                                                                                                |
| Interface             | 2x USB Type A<br>1x USB Type C                                                                                                                                                                                                                                                                                                                                                                                           |
| Weight                | 6.5 kg (14.3 lbs)                                                                                                                                                                                                                                                                                                                                                                                                        |
| Power supply          | 24 V  , max. 1.25 A                                                                                                                                                                                                                                                                                                                   |
| External power supply | Input: 100 - 240 V  ; 50 - 60 Hz ; max. 0.7 A<br>Output: 24 V  ; max. 1.25 A<br>Vendor: Mean Well<br>Model: GE30I24<br>Safety standards: UL62368-1, CSA C22.2 NO. 62368-1, TUV BS EN/EN62368-1, AS/NZS 60950.1, CCC GB4943, EAC TP TC 004 approved |

Technical data are subject to change without notice.

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**BYK-Gardner GmbH**